

Drug Testing Kits from Agriyork400 Ltd


Choose from Oral Fluid and Urine Multi Panel Drug Tests

Why Test?

Drug and alcohol use and abuse affects all of society. By monitoring, effective action and assistance can be provided.

A test programme acts as a deterrent.

Without monitoring, nothing is known until an incident has occurred.



Incidents caused by drug and alcohol impairment are NOT ACCIDENTS!

Correct testing combined with a modern and relevant drug panel make an effective Drug and Alcohol test programme possible.

Analyte	Result	Control
AMP	negative	23.05
COC	negative	30.32
MET	negative	34.81
MOR	negative	6.95
THC	negative	25.28

Why Oral Fluid?

Simple to sample – Less than 3 minutes.

No hygiene issues.

Short detection times = Recent Use.

Better indication of current state.

Observed sample collection means adulteration impossible.

Why Urine?

Longer detection times - 1 to 3 days

Adequate sample volumes for further confirmatory testing.

Why Hair?

In some cases, complete abstinence from illicit drugs and alcohol is demanded.

Hair analysis offers a long-term profile for drug detection.

CE Marked Drug tests available include:

Amphetamine, Barbiturates, Benzodiazepines, Buprenorphine, Cocaine, THC, Methadone, Methamphetamine, MDMA, Opiates, Oxycodone, PCP, Propoxyphene, TCA, Tramadol, EDDP, Ketamine, Fentanyl, K2/Spice

Which test panel is best?

A multi panel test makes more economic sense.

Drug tests with multiple test strips offer better detection rates.

Our AgriSAL Oral Fluid test is an 8 drug plus Alcohol Device (Picture below of a negative test result).

Our Euro200 Urine Cup test is a 12 panel Device.

The individual drugs on the test device you buy should ideally reflect the actual drug use and abuse in your community.

There is no point in paying and testing for drugs which are no longer used by people who take drugs in your area.

Testing for drugs which no longer reflect current AND local drug usage is a waste of time and resources.

How does it work?

The chemistry of the test is known as lateral flow immunoassay.

The sample, (urine or saliva) flows along a test strip.

The strip contains specially manufactured regions that will result in specific drugs reacting in test lines within the strip.

How do I read the results?

When the test is triggered by the flow of liquid (either urine or oral fluid) a Test T and Control C line are produced.

Control Line C forming means the test has successfully taken place.

Each drug has its own Test and Control Line.

Picture of a negative AgriSAL 8 panel test adjacent.

Easy to use. Training and certification available.

Reader Option Available for result storage and retrieval.

We monitor drug trends resulting in a current panel.

Results in minutes 5-minute development time.

Locally stocked warehouse, ensuring prompt shipments.

Long Shelf Life Typically 20 months.



Contact Agriyork400 Ltd at 01759 306580 or email sales@agriyork.co.uk
View our dedicated workplace site at www.workplacedrug.co.uk